

DNA IQ™ Method of Extracting DNA from Casework and Reference Samples

1. PURPOSE AND SCOPE	2
2. DEFINITIONS	2
3. PRINCIPLE	2
4. REAGENTS AND EQUIPMENT	4
4.1. Reagents	4
4.2. Extraction Buffer	5
4.3. Lysis Buffer with DTT Solution	6
4.4. DNA IQ™ Resin.....	6
4.5. 1x Wash Buffer	6
4.6. Elution Buffer	6
5. Equipment	7
6. SAFETY	7
7. SAMPLING AND SAMPLE PREPARATION	8
7.1. Sample Locations	8
7.2. QC Samples	8
7.2.1. Registration of QC Samples	8
7.3. Create the DNA IQ™ Lysis or Retain Supernatant batch	8
7.4. Locating Samples	8
8. OFF-DECK LYSIS PROCEDURE	8
8.1. Off-Deck Lysis (No Retained Supernatant)	8
8.2. Off-Deck Lysis (Retained Supernatant).....	10
9. AUTOMATED EXTRACTION OF LYSED SAMPLES	11
9.1. Create the DNA IQ Extraction batch	11
9.2. Locating samples	11
9.3. Sequence checking the Nunc Bank-It™ tubes	11
9.4. MPII Extraction Procedure	11
9.5. Summary of DNA IQ™ Extraction Version 6.5_ODL (following off-deck lysis)	11
9.6. Preparation of reagents for the automated extraction process	12
9.7. Setting up the MPII platforms for automated DNA IQ™ processing	13
9.8. Finalising the MP II Run	16
9.9. Importing MP II Log File into AUSLAB	16
9.10. Importing Extraction “results” into AUSLAB	17
9.11. Sample Storage	17
10. TROUBLESHOOTING	17
11. VALIDATION	17
12. QUALITY ASSURANCE/ACCEPTANCE CRITERIA	18
13. REFERENCES	18
14. STORAGE OF DOCUMENTS	19
15. ASSOCIATED DOCUMENTS	19
16. AMENDMENT HISTORY	19
17. APPENDIX	21
17.1. Manual method for extraction using DNA IQ™	21
17.1.1. Sampling and Sample Preparation	21
17.1.2. QC samples	21
17.1.3. Creating the Extraction Batch and Locating Samples	21

Automated DNA IQ™ Method of Extracting DNA

17.1.4. Procedure (No Retain Supernatant).....	21
17.1.5. Procedure (Retain Supernatant).....	23
17.1.6. Sample storage.....	25

1. PURPOSE AND SCOPE

This method describes the routine method for the extraction of DNA using the DNA IQ™ kit (Promega Corp., Madison, WI, USA). The automated method is the preferred procedure, utilising the MultiPROBE® II PLUS HT EX with Gripper™ Integration platforms (PerkinElmer BioDiscovery, Downers Grove, IL, USA). The manual method has also been included. This method applies to all DNA Analysis staff members that are required to extract DNA from samples.

Reference samples and casework samples must be extracted separately. If casework and reference samples are to be extracted on the same instrument, the instrument (including all required labware) must be decontaminated between operations.

2. DEFINITIONS

DNA IQ™ Resin	Magnetic resin beads used to bind DNA
DTT	1,4 Dithiothreitol
EDTA	Ethylenediaminetetraacetic acid
EP-A	Extraction Platform A
EP-B	Extraction Platform B
Extracts	Samples that have had a DNA extraction processes performed
Lysates	Samples that have had the off-deck lysis step performed, but have not yet completed the entire extraction process
MPII	MultiPROBE® II PLUS HT EX Platform
Paramagnetic	To become magnetic with the application of a magnetic force
Pro K	Proteinase K
Samples	Sample substrates (in tubes) awaiting DNA extraction
Sarcosyl	N-Lauroylsarcosine sodium
TNE	Tris, NaCl and EDTA buffer (10mM Tris, 100mM NaCl, 1mM EDTA, pH 8.0)

3. PRINCIPLE**Sample Pre-lysis**

The Extraction Buffer used in the pre-lysis of cellular material contains Proteinase K (Pro K), TNE buffer (10mM Tris, 100mM NaCl, 1mM EDTA, pH 8.0) and Sarcosyl. TNE acts as a basic buffer with EDTA chelating ions in solution. Sarcosyl is a detergent that lyses open cell membranes. Proteinase K is added to digest cellular material, helping the DNA to be released. It also digests proteins that interfere with the DNA binding capacity of the resin. In addition, Proteinase K rapidly inactivates enzymatic activity that could potentially degrade DNA (e.g. nucleases).

Proteinase K (or Endoproteinase K) is a fungal enzyme from *Engyodontium album* (formerly *Tritirachium album*). It is a stable S8 family serine alkaline protease that has broad specificity and is inactivated by diisopropylfluorophosphate, phenyl-methane sulfonyl fluoride and Hg²⁺ ions. As a highly active endopeptidase, it cleaves peptide bonds adjacent to the carboxyl group of N-substituted (blocked alpha amino) hydrophobic, aliphatic and aromatic amino acids. Therefore, it actively hydrolyses proteins, glycoproteins, peptides and esters of amino acids. The smallest peptide that can be hydrolysed with Proteinase K is a tetrapeptide.

Automated DNA IQ™ Method of Extracting DNA

DNA IQ™ Kit

The DNA IQ™ kit (Promega Corp., Madison, WI, USA) is designed to purify DNA from a number of different substrates and has been optimised for trace DNA samples. An in-house validation was performed using a modified version of the PerkinElmer automated protocol.

The in-house protocol includes:

- Off-deck lysis steps with the option to retain a portion of the supernatant for further testing;
- The use of 300µL Extraction Buffer containing TNE (10mM Tris, 100mM NaCl, 1mM EDTA, pH 8.0) and Proteinase K to lyse cellular material prior to performing the DNA IQ process;
- The use of tubes and spin baskets for off-deck lysis of samples prior to extraction on the MPII platform. At the conclusion of off-deck lysis, lysates are transferred to individual Nunc Bank-It™ tubes;
- Nunc Bank-It™ tubes (arranged in sequence using STORstar) containing lysates are presented to the MPII platform for automated transfer of lysates into a 96-deep well plate;
- DNA IQ™ Resin is added using the MPII platform, followed by addition of two volumes of DNA IQ™ Lysis Buffer;
- The 96-deep well plate containing DNA IQ™ Resin and Lysis Buffer is sealed using an adhesive aluminium film and is placed on a MixMate to mix the contents of each well. The plate is centrifuged and the aluminium film is then pierced using a 96 well half skirt PCR microplate and the plate is returned to the MPII platform;
- A double elution step is performed using two dispenses of DNA IQ™ Elution Buffer at 60µL, resulting in a final DNA extract volume of 100µL;
- DNA extracts are automatically transferred into Nunc Bank-It™ tubes for storage.

Cell lysis is performed using DNA IQ™ Lysis Buffer containing Dithiothreitol (DTT). 1,4 Dithiothreitol is a reducing agent used in extractions to break disulfide bonds of proteins. The Lysis Buffer used is a proprietary buffer containing chaotropic salts required for binding of DNA to the magnetic beads. According to the MSDS (Promega Corporation 2006), the Lysis Buffer contains 50-75% guanidinium thiocyanate, < 2.5% EDTA, < 2% 3-[(3-Choalamidopropyl)dimethylammonio] propanesulfonic and < 2% polyethylene glycol tert-octylphenyl ether.

The basis of the DNA IQ™ kit is a magnetic bead resin which contains novel paramagnetic particles covered with silica. The magnetic bead resin usually has a DNA binding capacity of 100ng but the addition of Proteinase K increases the binding capacity. Samples with small amounts of DNA are more efficient than samples with large amounts of DNA, and as a result a small sample size is critical to ensure efficient recovery of DNA.

The magnetic beads with silica have a negative charge at basic and near neutral pH. The Lysis Buffer changes the pH and salt concentration of the solution and the silica on the magnetic beads becomes positively charged, which is then able to bind the DNA.

Several washing steps are employed in the protocol to remove inhibitors. The first washing step is performed using Lysis Buffer. This wash ensures the DNA is bound to the magnetic bead resin and washes out inhibitors. The next three washing procedures incorporate the use of DNA IQ™ Wash Buffer. This buffer contains an alcohol/aqueous mixture which ensures the DNA is not eluted during washing by keeping the DNA dehydrated, and the aqueous phase washes out the inhibitor.

The DNA IQ™ Elution Buffer removes the DNA from the magnetic bead resin. The Elution Buffer changes the salt content. Heating the complex to 65°C allows the DNA to be released from the magnetic beads. The Elution Buffer is a low ionic strength buffer that

Automated DNA IQ™ Method of Extracting DNA

reduces the affinity of the DNA for the silica that covers the magnetic beads, by re-hydration of the phosphate backbone.

The DNA IQ™ kit isolates DNA greater than 80bp, smaller DNA is removed selectively to prevent PCR inhibition.

MultiPROBE® II HT EX PLUS with Gripper™ Integration Platform

Within DNA Analysis, routine DNA extractions are performed on casework or reference samples using two MultiPROBE® II PLUS HT EX with Gripper™ Integration platforms (EP-A or EP-B) located in Room 6127.

Each platform uses a computer-controlled Cartesian XYZ liquid handling system designed for the efficient automation of sample preparation. Liquid transfer is performed by an 8-tip system with VersaTip® and VariSpan™ options. The VersaTip® option allows the use of both fixed and disposable tips (conductive and non-conductive). The VariSpan™ option permits variable probe spacing between each of the individual probes so that a wide range of labware such as micro plates, tubes, vials and reagent troughs can be accessed. Each sample probe is capable of independent motion in the Z direction due to independent Z drives.

The 8-tip system is also capable of multichannel liquid-level sensing by utilising Accusense™ technology. This technology works by each probe detecting a change in capacitance within the liquid. This capacitive mode of detection is also possible when using conductive disposable tips. Accusense™ also permits the detection of low or non-ionic polar solutions and solvents. Pipetting on the platforms is driven by 8 individually controlled syringes via the positive displacement of the system liquid (nanopure water) when transferring liquid.

The Gripper™ Integration allows for automated identification of labware via the scanning of barcodes and also allows for automated transfer of plates from one position on the deck to another. To increase deck space and capacity, the platforms include a left deck extension.

For automated DNA extraction using the DNA IQ™ kit, a platemap is utilised to provide the necessary information for correct volumes and locations for pipetting. It also contains information regarding batch and plate identification, sample positions, lab numbers of samples, and lists all the samples in the batch, providing the total number of samples. The program will fail to work without a platemap.

4. REAGENTS AND EQUIPMENT

4.1. Reagents

- DNA IQ™ System Kit (400 sample kit)
 - DNA IQ™ Resin
 - Lysis Buffer (LB)
 - 2x Wash Buffer (2xWB)
 - Elution Buffer (EB)
- TNE (10mM Tris, 100mM NaCl, 1mM EDTA, pH 8.0)
- Proteinase K (Pro K) 20mg/mL
- Dithiothreitol (DTT) 1M
- 5% TriGene
- 70% Ethanol
- 10% Bleach 7x Solution
- 1% Amphyl
- 0.2% Amphyl
- Isopropyl Alcohol
- AnalR 100% Ethanol

Automated DNA IQ™ Method of Extracting DNA

- 40% Sarcosyl
- Nanopure Water

These reagents are stored in locations as per Table 1.

Table 1. Reagent storage locations.

Reagent	Device	Storage Location
Pro K	Freezer	Room 6120
DTT	Freezer	Room 6120
40% Sarcosyl	Shelf	Room 6122
Isopropyl Alcohol	Shelf	Room 6122
AnalR 100 %Ethanol	Shelf	Room 6122
TNE pH 8 Buffer	Shelf	Room 6122
DNA IQ™ Kit	Shelf	Room 6122
Amphyl (1% and 0.2%)	Shelf	Room 6127
Nanopure Water	Shelf	Room 6127
5% TriGene	Shelf	Room 6127
10% Bleach 7x Solution	Shelf	Room 6127

Table 2 shows the volume of reagents for a full plate or half plate. All reagents can be made on the bench, with the exception of the Lysis Buffer-DTT which needs to be made in a fume hood. DNA IQ™ reagents are prepared by staff performing the method. Refer to “*Receipt, Storage and Preparation of Chemicals, Reagents and Kits*” (QIS [17165](#)) for preparation of TNE Buffer.

Table 2. Table of reagent volumes.

Reagent	Volume for 96 samples (mL)	Volume for 48 samples (mL)	Volume for 24 samples (mL)
Extraction Buffer			
TNE buffer	33.3	20	10
Proteinase K (20mg/mL)	1.8	1.08	0.54
Sarcosyl (40%)	0.9	0.54	0.27
Lysis-DTT Buffer			
DNA IQ™ Lysis Buffer	90.0	50	N/A
DTT (1M)	0.9	0.5	N/A
DNA IQ™ Resin solution			
Lysis-DTT Buffer	6.0	3	N/A
DNA IQ™ Resin	1.0	0.5	N/A
DNA IQ™ 1x Wash Buffer			
DNA IQ™ 1x Wash Buffer	35.0	18	N/A
DNA IQ™ Elution Buffer			
DNA IQ™ Elution Buffer	14.0	8	N/A

Note: Do not directly pipette from reagent stock bottles. Stock reagents should be decanted into working containers where appropriate. Volume for 24 samples is for off-deck lysis samples only.

4.2. Extraction Buffer

Note: Prepare Extraction Buffer just prior to commencing the off-deck lysis or extraction procedure.

1. Determine the required volumes of reagents by using Table 2.
2. Remove the required amount of 20mg/mL Proteinase K from the freezer and thaw. Vortex and centrifuge before use.
3. Ensure that the 40% (w/v) Sarcosyl is completely dissolved and appears clear in the stock solution. If not dissolved, invert the container a few times and leave at room temperature.
4. Retrieve an aliquot of TNE buffer of the appropriate volume size from the falcon tube storage container in Room 6122.
5. Add the appropriate volumes of 20mg/mL Proteinase K and 40% (w/v) Sarcosyl to the falcon tube containing TNE buffer, and invert gently to mix.
6. Label the tube with “Extraction Buffer”, your initials and the date.

4.3. Lysis Buffer with DTT Solution

Note: Lysis Buffer is supplied with the DNA IQ™ kit. The Lysis Buffer with DTT solution is prepared just prior to commencing the extraction procedure.

Warning: Lysis Buffer and DTT are toxic, use a fume hood to prepare these reagents and wear appropriate PPE including safety glasses when handling Lysis Buffer.

1. Determine whether a half- or full-plate of reagents are required (Table 2).
2. Remove the required amount of DTT from the freezer and thaw. Vortex and centrifuge before use.
3. In the fume hood add the required volume of Lysis Buffer to a sterilised glass Schott bottle and then add the required volume of DTT.
4. Label the glass Schott bottle with "Lysis Buffer + DTT", your initials and the date.

4.4. DNA IQ™ Resin

Note: DNA IQ™ Resin is supplied with the DNA IQ™ kit. The resin is prepared at the start of each run. Ensure the resin is properly mixed by *vortexing* prior to use.

1. Determine whether a half- or full-plate of reagents are required (Table 2).
2. Into a 10mL (or 5mL) sterile tube, add the required volume of Lysis Buffer with DTT solution (from 4.1.2) followed by the required volume of DNA IQ™ Resin.
3. Mix by gentle inversion.
4. Label the tube with "Resin", your initials and the date.

4.5. 1x Wash Buffer

Note: 2x Wash Buffer is supplied with the DNA IQ™ kit. To prepare 1x Wash Buffer, add 35mL of *AnalR* Ethanol and 35mL of Isopropyl Alcohol to the 2x Wash Buffer bottle. Then label the lid and side of the bottle with "1x Wash Buffer," your initials and the date. Also fill out the Reagent Log (DNA IQ Reagents).

1. Determine whether a half- or full-plate of reagents are required (Table 2).
2. Into a Falcon tube, add the required volume of 1x Wash Buffer.
3. Label the falcon tube with "Wash Buffer", your initials and the date.

4.6. Elution Buffer

Note: Elution Buffer is supplied with the DNA IQ™ kit. The Elution Buffer can be used directly from the kit. The Elution Buffer is removed from the kit and stored in the automated extraction room (6127).

Automated DNA IQ™ Method of Extracting DNA

5. EQUIPMENT

The following equipment (Table 3) and consumables (Table 4) are required for the DNA IQ™ extraction process.

Table 3. Equipment used and location.

Equipment	Asset No	Location
STORstar (B system)		6122
MultiPROBE® II PLUS HT EX with Gripper™ Integration Platform (EP-A)		6127
MultiPROBE® II PLUS HT EX with Gripper™ Integration Platform (EP-B)		6127
DPC shaker (EP-A)		6127
DPC shaker (EP-B)		6127
Automated Temperature Controller, Heat Block tiles and heat block adapters (EP-A)		6127
Automated Temperature Controller, Heat Block tiles and heat block adapters (EP-B)		6127
Eppendorf 5804 centrifuge		6127
Vortex		6127
Fridge		6127
Micro centrifuge		6127
MixMate		6127
Decapper		6127
4itude 4seal Sealer		6127

Table 4. Consumables used for extraction.

Consumables	Location
175µL Clear Non-Conductive Filter RoboRack tips – Pre-Sterilised	6127
MβP Pure 1000uL Tips – Pre-Sterilised	6127
SlicPrep™ 96 device plate	6122
ABgene 96-deep well plate	6120
Axygen 2mL deep well storage plate	6127
96 well Half Skirt PCR Microplate	6127
1.5mL or 2mL Eppendorf tubes with spin baskets	6120
12 Channel plate	6127
Nunc Bank-it™ tubes	6120
Nunc Bank-it™ caps	6127
Sterile 50mL Falcon tubes	6122
Sterile 10mL or 5mL tubes	6122
Autoclaved 100mL glass bottles	6122
Autoclaved 250mL glass bottles	6122
Aluminium sealing film	6127
300µL ART tips	6120
1000µL ART tips	6120

6. SAFETY

As per the procedures in the QIS document “Operational Practices in the DNA Dedicated Laboratories” (QIS 17120), PPE is to be worn by all staff when performing this procedure. This includes the use of safety glasses.

The automated platforms, labware and associated equipment should be wiped with 5% TriGene™ followed by 70% Ethanol before and after use. In addition, clean the work area around the MPII instrument with 10% (v/v) Bleach 7x and 70% Ethanol.

While the MPII is running a procedure, refrain from placing any part of the body inside the cabinet. Either pause the program or use the emergency STOP button located at the front of the instrument before touching anything on the deck surface. Pressing the emergency STOP button may cause the program to pause or abort.

Warning: Tris base, EDTA, Sarcosyl and Lysis Buffer are irritants. DTT is a reducing agent that destroys disulfide bonds. Handle carefully and wear appropriate PPE.

Lysis Buffer contains guanidinium thiocyanate (GuSCN) which is toxic and can be harmful if inhaled, swallowed or comes in contact with skin. Left over Lysis Buffer-DTT is disposed of

Automated DNA IQ™ Method of Extracting DNA

in a brown Winchester bottle in the fume hood. Never dispose down the sink. If spillage occurs, absorb with liquid-binding material (such as sand or kitty litter) and dispose in a biohazard bin. Do not combine bleach with Lysis Buffer and do not wipe areas of Lysis Buffer spillage with bleach, as the two reagents produce cyanide gas when combined. Handle carefully and wear appropriate PPE, including safety glasses, when handling Lysis Buffer. If Lysis Buffer is spilt onto PPE (eg. gloves, gowns), discard the PPE and obtain new PPE.

7. SAMPLING AND SAMPLE PREPARATION

7.1. Sample Locations

Samples waiting to be extracted are stored in freezers as described in Table 5.

Table 5. Sample storage locations.

Sample type	Storage Device	Storage Location
Urgent/High/Medium Priority Samples	Freezer	6117-2
Lysates in 1.0mL Nunc Bank-it™ tubes	Freezer	6117
Extracts in 1.0mL Nunc Bank-It™ tubes	Freezer	6117

Note: Some Medium and Low Priority storage boxes may be located in the Exhibit Room (6106).

7.2. QC Samples

For all off-deck lysis batches (with 48 samples or less) and extraction batches; one negative control and one positive control is required to be registered. For all off-deck lysis batches with > 48 samples; two negative and two positive controls is required to be registered.

Table 6. Extraction Quality Controls

QC Name	Batch Type	Description
Negative Control	Off-Deck Lysis	Negative Extraction control – empty well
Positive Control	Off-Deck Lysis	Positive Extraction control – dried blood swab from a known donor
Negative Control	IQ Extraction	Negative Extraction control – empty well
Positive Control	IQ Extraction	Internal IQ Efficiency Control

7.2.1. Registration of QC Samples

The registration of control samples is covered in the DNA Analysis workflow procedure (QIS [24919](#))

7.3. Create the DNA IQ™ Lysis or Retain Supernatant batch

Creation of Lysis and retain supernatant batches is covered in the DNA Analysis Workflow Procedure (QIS 24919).

7.4. Locating Samples

To locate samples refer to “Analytical Sample Storage” (QIS [24255](#)).

8. OFF-DECK LYSIS PROCEDURE

8.1. Off-Deck Lysis (No Retained Supernatant)

- For batches of 48 samples or less, one set of controls and one operator is required. For larger batches, separate the batch into two smaller batches of 48 samples, including one set of controls in each. If a single operator is performing the whole procedure, the second batch can be started during an appropriate incubation step.

Note: For full batches positions 1-4 will be the two sets of controls, positions 5-50 will be the first 46 samples and 51-96 will be the second set of 46 samples.

Automated DNA IQ™ Method of Extracting DNA

2. For each sample label:
 - Original sample tube
 - Spin basket or 1.5mL tube as required
 - 1.0mL Nunc Bank-It™ tube

Note 1: Samples that require a spin basket are: swabs, fabric, paper and cigarette butts. Samples that require a 1.5mL tube are: tapelifts, chewing gum, straws, fingernail clippings/scrapings and toothbrush bristles.

Note 2: If samples are in a 2mL QPS tube and require a spin basket, label a new tube for the substrate to be retained in.

3. If substrate is in a 0.5mL tube, transfer to a labelled 1.5mL tube.
4. Have a second operator sequence check the tubes and enter the sequence check in AUSLAB.
5. Prepare Extraction Buffer as per Section 4.1.1.
6. Add 300µL of Extraction Buffer and vortex briefly. Ensure that substrates are immersed in the Extraction Buffer.
7. Incubate on a Thermomixer at 37°C for 45 minutes at 1000rpm or on the hotblock (if using hotblock vortex samples at 22-23 minutes and again at the end of the incubation). Record temperature on worksheet.
8. Remove from the Thermomixer/hotblock. Transfer substrate to spin basket if required or transfer the lysate to an appropriately labelled 1.5mL tube. Retain original tube containing the substrate in if no spin basket used.
9. Centrifuge spin basket at maximum speed (14,000rpm) for 2 minutes.
10. Retain spin basket containing the substrate and transfer flow through back to original lysis tube.
11. Vortex lysate, then incubate in hotblock/Thermomixer at 65°C for 10 minutes. Record temperature on worksheet.
12. Centrifuge at maximum speed (14,000rpm) for 1 minute.
13. Transfer 300uL of lysate to the corresponding Nunc Bank-It™ tube.

Note: If more than 300uL of lysate is present, retain the remaining lysate in the 1.5mL tube. In AUSLAB, write a specimen note for that sample (eg. "extra lysate retained from sample XXXXXXXXX"). Store the retained 1.5mL lysate tube in appropriate box in freezer.

14. In AUSLAB, enter reagent and hotblock/Thermomixer temperature details and complete the batch.
15. Transfer substrates from spin baskets to an appropriately labelled tube (may use original sample tube if no remaining lysate)
16. Store lysates in temporary storage boxes in freezer 6117-2 (-20°C). Store tubes containing substrates in "Spin Basket boxes" in freezer 6117-5 (-20°C).

*Automated DNA IQ™ Method of Extracting DNA***8.2. Off-Deck Lysis (Retained Supernatant)**

1. For batches of 48 samples or less, one set of controls and one operator is required. For larger batches, separate the batch into two smaller batches of 48 samples, including one set of controls in each. If a single operator is performing the whole procedure, the second batch can be started during step 7.

Note: For full batches positions 1-4 will be the two sets of controls, positions 5-50 will be the first 46 samples and 51-96 will be the second set of 46 samples.

2. For each sample label:
 - Original sample tube
 - Spin basket or 1.5mL tube as required
 - 1.5mL tube (also labelled with “sup” to indicate supernatant)
 - 1.0mL Nunc Bank-It™ tube

Note 1: Samples that require a spin basket are: swabs, fabric, paper and cigarette butts. Samples that require a 1.5mL tube are: tapelifts, chewing gum, straws, fingernail clippings/scrapings and toothbrush bristles.

Note 2: If samples are in a 2mL QPS tube and require a spin basket, label a new 5mL tube for the substrate to be retained in.

3. If substrate is in a 0.5mL tube, transfer to a labelled 1.5mL tube.
4. Have a second operator sequence check the tubes and enter the sequence check in AUSLAB.
5. Add 450µL of TNE buffer and vortex.
6. Incubate at room temperature for 30 minutes.
7. Vortex, then centrifuge at maximum speed (14,000rpm) for 3 minutes.
8. Remove 150µL of supernatant and place into the respective 1.5mL tube labelled with “sup” (for further testing).
9. Add 14µL of 20mg/mL Proteinase K and 7µL 40% (w/v) Sarcosyl to the original sample tube containing substrate and TNE buffer, then vortex.
10. Incubate in Thermomixer at 37°C for 45 minutes at 1000rpm or on the hotblock (if using hotblock vortex samples at 22-23 minutes and again at the end of the incubation). Record temperature on worksheet.
11. Remove from the Thermomixer/hotblock. Transfer substrate to spin basket if required or transfer the lysate to an appropriately labelled 1.5mL tube. Retain original tube containing the substrate in if no spin basket used.
12. Centrifuge spin basket at maximum speed (14,000rpm) for 2 minutes.
13. Retain spin basket containing the substrate and transfer flow through back to original lysis tube.
14. Vortex Lysate, then incubate in hotblock/Thermomixer at 65°C for 10 minutes. Record temperature on worksheet.
15. Centrifuge at maximum speed (14,000rpm) for 1 minute.

Automated DNA IQ™ Method of Extracting DNA

16. Transfer 300uL of lysate to the corresponding Nunc Bank-It™ tube.

Note: If more than 300uL of lysate is present, retain the remaining lysate in the 1.5mL tube. In AUSLAB, write a specimen note for that sample (eg. “extra lysate retained from sample XXXXXXXXX.”). Store the retained 1.5mL lysate tube in appropriate box in freezer.

17. In AUSLAB, enter reagent and hotblock/Thermomixer temperature details and complete the batch.

18. Transfer substrates from spin baskets to an appropriately labelled tube (may use original sample tube if no remaining lysate)

19. Store supernatants in the “S/N Retention” boxes in Freezer 6117-2 (-20°C). Store lysates in temporary storage boxes in freezer 6117-2 (-20°C). Store tubes containing substrates in “Spin Basket boxes” in freezer 6117-5 (-20°C).

9. AUTOMATED EXTRACTION OF LYSED SAMPLES

9.1. Create the DNA IQ Extraction batch

Creation of extraction batch is covered in the DNA Analysis Workflow Procedure (QIS 24919).

9.2. Locating samples

To locate samples refer to “Analytical Sample Storage” (QIS [24255](#)).

9.3. Sequence checking the Nunc Bank-It™ tubes

The procedure for the automated checking of sample tubes is covered in the Procedure for the use of the STORstar unit for automated sequence checking (QIS 24256)

9.4. MPII Extraction Procedure

Automated setup of DNA IQ extractions in 96-deep well format is performed using the MultiPROBE® II PLUS HT EX EP-A and EP-B platforms located in Room 6127.

Refer to “*Operation and Maintenance of the MultiPROBE® II PLUS HT EX and MultiPROBE® II PLUS HT EX with Gripper™ Integration Platform*” (QIS [23939](#)) for instructions on the use and maintenance of the MultiPROBE® II PLUS HT EX platforms.

9.5. Summary of DNA IQ™ Extraction Version 6.5_ODL (following off-deck lysis)

1. Transfer of lysates from Nunc Bank-It™ tubes into the ABgene 96-deep well plate

Lysates from the off-deck lysis protocol, contained in individual Nunc Bank-It™ tubes, are transferred automatically into an ABgene 96-deep well plate prior to commencing automated sample processing. Instead of a 96-deep well plate, the use of individual Bank-It™ tubes for storage of lysates removes the need for an adhesive film and therefore reduces the risk of well-to-well contamination from the action of peeling back an adhesive plate cover.

2. Automated addition of DNA IQ™ Resin and Lysis Buffer

DNA IQ™ Resin is added automatically into the ABgene 96-deep well plate on the platform. The program uses a waste mode dispensing step to add 53µL DNA IQ™ Resin-Lysis Buffer solution to each well in either half the plate or the whole plate. Because of this, batch sizes are restricted to either 48 or 96 samples in order to

Automated DNA IQ™ Method of Extracting DNA

maintain efficiency and economy. Two volumes of Lysis Buffer are then added to promote binding of DNA to the paramagnetic resin.

3. Mixing using a MixMate to bind DNA to resin

Manual intervention is required to seal the ABgene 96-deep well plate with a 4titude Pierce Seal and sealing plate using 4titude sealer pre-heated at 175 °C. The plate is then transferred onto a MixMate instrument for mixing at 1100rpm for 5 minutes, followed by centrifugation on an Eppendorf 5804 centrifuge at 3000rpm for 2 minutes. The seal is then carefully pierced with a 96 well PCR microplate and the ABgene plate is returned to the Applied Biosystems magnet on the MPII platform.

4. Removing lysis reagents for storage

At this point, most of the DNA is bound to the paramagnetic resin. With the positioning of the ABgene plate on the ABI magnet, DNA IQ™ Resin becomes immobile at the bottom of the plate. The lysis reagents from each well are transferred automatically to a storage plate on the MPII platform without disturbing the DNA IQ™ Resin. The purpose of the storage plate is for retaining supernatant that may potentially still contain DNA material. The storage plate may also become useful in quality investigations.

5. Washing of the resin-DNA complex

Washing steps are performed to remove any inhibitors in solution. The first wash uses 125µL Lysis Buffer with shaking at room temperature for 1 minute on the DPC shaker to ensure that the DNA is bound to the paramagnetic resin. The plate is moved to the ABI magnet and the supernatant is transferred into the storage plate. The next three washes incorporate 100µL of 1x Wash Buffer with shaking at room temperature for 1 minute on the DPC shaker. During each wash cycle, the ABgene 96-deep well plate is moved to the ABI magnet and the supernatant is discarded into the tip chute.

6. Removing any excess of 1x Wash Buffer

The samples are allowed to air dry at room temperature for 5 minutes in order to remove all traces of Ethanol from the Wash Buffer. The presence of Ethanol may potentially inhibit both the elution process and also downstream PCR.

7. Elution of DNA from the resin-DNA complex

A double elution method is employed in this procedure. At each elution step, 60µL of Elution Buffer is added to each sample, followed by incubation at 65°C for 6 minutes (3 minutes without shaking followed by 3 minutes shaking on the DPC shaker). The ABgene plate is moved to the ABI magnet and the eluted solution (supernatant containing eluted DNA) is transferred to fresh Nunc™ Bank-It™ tubes.

8. Flushing of capillaries

As a decontamination measure, the MPII capillaries and liquid pathway are washed with Amphyl and Nanopure water at the conclusion of the automated extraction process.

9.6. Preparation of reagents for the automated extraction process

Note: Reagents are prepared during the setting up of the MPII platforms (Section 4.3).

9.7. Setting up the MPII platforms for automated DNA IQ™ processing


The following steps are carried out in the automated extraction room (Room 6127).

1. Remove the Nunc Bank-It™ tubes containing lysates from the fridge to allow to come to room temperature before commencing the extraction procedure.

Note: If the lysates are frozen, remove them from the freezer and thaw in Room 6127. Also remove the required amount of DTT to thaw.

2. Restart or turn on the instrument PC.
3. Log onto the network using the **Robotics** login.
4. Open WinPrep® by double clicking icon on the computer desktop (Figure 1).
5. Log onto the WinPrep® software by entering your username and password, then press “Enter”.
6. Ensure that the daily/weekly start-up and maintenance has been performed before running any program. If WinPrep® has been closed or been idle for a long period of time initialise the MP II platform as described in QIS 23939.
7. Ensure the **System Liquid reservoir is FULL** and tubing is fully submerged in the system liquid before every run and perform a Flush/Wash. If visible air bubbles have appeared in tubing or in syringes between setting up the deck and executing the test, another flush wash will need to be performed before starting the MPII extraction run.
8. Open the Extraction setup MP II test file in WinPrep® by selecting:
 - **File**
 - **Open**, navigate to **C:\PACKARD\MULTIPROBE\BIN\QHSS PROTOCOLS**
 - Select “**DNA IQ Extraction_Ver 6.5_ODL.mpt**”
 - Click the “**Open**” button
9. Check the tree pane of the whole program for any bold fonts. See the Analytical Senior Scientist if bold fonts are present.
10. Copy the required plate map from the network folder **I:\EXTRACTION** into the local folder **C:\PACKARD\EXT PLATE MAPS**. Open the plate map using Excel and check that the plate map is complete, and make changes if necessary. If there are blanks (samples that have been removed), make all volumes to 0. If it is not a full batch, delete the rows of blank samples at the end of the platemap. Save all changes made to the platemap (as a .txt file).
11. Decontaminate the required labware with 5% TriGene followed by 70% Ethanol and place onto the instrument deck in the designated grid positions (as shown in the WinPrep® software). Additionally, ensure the DPC shaker is positioned properly.
12. Ensure that the DPC shaker and Heater Controller Box are switched on.
 - For EP-A: Tile 1 at F22 (85°C).
 - For EP-B: Tile 2 at F22 (85°C).

Note: Press the start/stop button twice at the front of the DPC shaker to ensure that it displays zero on the screen.



Automated DNA IQ™ Method of Extracting DNA

13. Ensure the heat transfer tile is clicked into the plate adapter tile properly.
Note: This is critical to ensure correct incubation temperatures.
14. To the Amphyl wash station in position **A10**, add 1% Amphyl to the far left hand side, 0.2% diluted Amphyl to the middle left and place an empty (waste) reagent trough in the middle right position of the station.
15. Refer to section [4.1](#) for the preparation of reagents. Record all lot numbers onto the worksheet and in AUSLAB. Note, for batches of <48 samples, use volumes for 48 samples.
16. Check the syringes and tubing and perform a Flush/Wash if required.
17. Add Wash Buffer to the far right hand side trough of the Amphyl wash station in position **A10**. Add Lysis Buffer to the labelled 150mL reagent trough on the left hand side of the 2 trough holder in position **A13**.
18. Place the 12 channel plate into position **A16**. Add Elution Buffer to the plate by dividing the Elution Buffer between channels 11 and 12. Add Resin to channel 1. It is important to add the resin in a uniform fashion to ensure equal distribution of resin along the channel.
19. **Nunc Bank-It™ lysate tubes:** The lysates should now be at room temperature. Ensure that the rack is labelled with the correct **AUSLAB batch ID** on the **front** of the Nunc™ Bank-It™ tube rack and that the label matches the batch ID on the worksheet and platemap. Additionally, randomly check some visible barcodes on the side of lysate tubes against the positions on the worksheet to ensure the correct batch of samples is to be processed.
 - a. Add a B1-Lite generated '**LYSATE**' barcode on the **right hand side** of the Nunc™ Bank-It™ tube rack.
 - b. Place the rack onto the MixMate to shake for 1 minute at 1000rpm.
 - c. Centrifuge the rack at 3000rpm for 2 minutes in the Eppendorf 5804 centrifuge and then place into position **C13**.**Note:** Do not uncap lids until prompted by program.
20. **ABgene 96-deep well plate:** Label the **left hand side** of the plate with both the correct **AUSLAB batch ID** and **batch ID barcode**. With a marker, print the word "FRONT-lysate" on the front of the plate. Place the plate in its correct orientation in position **E13**.
21. **2mL 96-deep well storage plate:** Label the **left hand side** of the plate with both the correct **AUSLAB batch ID** and **batch ID barcode**. Label the **right hand side** of the plate with a B1-Lite generated "**STORE**" barcode. With a marker, print the word "FRONT-store" on the front of the plate. Place the plate in its correct orientation in position **E16**.
22. **Nunc Bank-It™ extract tubes:** Ensure that the rack is labelled with the correct **AUSLAB batch ID** on the **front** of the Nunc™ Bank-It™ tube rack. Label the **right hand side** of the plate with a B1-Lite generated "**EXTRACT**" barcode. Additionally, randomly check some visible barcodes on the side of lysate tubes against the positions on the worksheet to ensure the correct batch of samples is to be processed. Place the rack in position **G16**.
Note 1: Do not uncap lids during this step.
Note 2: If B1-Lite generated barcodes are not available hand-write the labels.
23. Add Nanopure water to the 160mL trough in the Flush/Wash station in position **G13**.

Automated DNA IQ™ Method of Extracting DNA

24. Ensure that all necessary labware have been positioned correctly as displayed within WinPrep®, then click **“EXECUTE TEST”**. Record run information in the Run Log book.
25. The following message will appear (Figure 2 below):




Figure 2. Scan batch ID request

Into “New Rack ID:” scan barcode off the worksheet. It is important this corresponds to the labelling of labware on the deck and the plate maps used.

26. Click **“Reset Tip Boxes”** and ensure that the number of tips displayed in this window match those as available on the deck platform. Fill new tips if necessary. Click **“Close”** to accept the tip count, and then click **“Next”**.
27. Select the correct platemap by browsing to **C:\PACKARD\EXT PLATE MAPS**. Ensure that the platemap selected corresponds to the labelling of labware on the deck, the paperwork used and the “New Rack ID” entered above.
28. For a full batch of 96 samples, ensure that all nodes are checked. For a batch of 48 samples or less, uncheck the node: “Add resin to second half of plate”. Click **“Next”** to check all other nodes.
29. Click **“Start”** to continue.
30. The MPII instrument will proceed to scan the required plates on the platform deck in the below order. If barcode reading fails or if B1-Lite barcodes are not available (and hand-written labels have been used), the user is prompted to enter a plate ID. A plate ID can be entered manually into the “Read failed” prompt window for:
- Nunc extract tubes, type in **EXTRACT** and press **“Enter”**.
 - 96-deep well storage plate, type in **STORE** and press **“Enter”**.
 - Nunc lysate tubes, type in **LYSATE** and press **“Enter”**.
31. After the plates have been identified, two user prompts will appear as a reminder to confirm the deck setup. **Always decap tubes from positions H1 to A1, H2 to A2 etc.**
- Ensure all steps on the first prompt have been complete, Click **OK** to continue.
 - Ensure all steps on the second prompt have been complete, Click **OK** to continue.
32. The program will progress to transfer the lysates followed by automated addition of the Lysis-DTT buffer and the DNA IQ™ Resin solution. The next user prompt will appear. Follow the steps as outlined in the user prompt. Then click **OK** when ready.
Note: Ensure that plate is sealed properly with the Pierce Seal. Once the Pierce Seal film is pierced, the PCR Microplate is then discarded (new plate used each time).
33. Once lysis steps have been completed, a user prompt will appear. Follow the directions as outlined in the user prompt. Then click **OK** when ready.
Note: Nunc lysate tubes (including the rack) should be placed into a clipseal bag, sealed and discarded into a biohazard waste bin.

Automated DNA IQ™ Method of Extracting DNA

34. Once the wash steps have been completed, a user prompt will appear. Follow the directions as outlined in the user prompt. Then click **OK** when ready.
Note: The Nunc extract tubes can be decapped on the bench in the MP11 cabinet while the samples are drying and the heating tile is stabilising. At this point, the operator has approximately **12 minutes** to decap the tubes. If the operator has been interrupted and requires more time to decap the tubes, the program may be paused at this step.
35. A message will appear waiting for the heating tile to reach 85°C (for incubation at 65°C). **DO NOT PRESS CONTINUE** as the program will continue automatically when the temperature has been reached with sufficient stability.
36. A user prompt will appear. Follow the directions as outlined in the user prompt. Then click OK when ready.
37. Once the elution is completed, a user prompt appears. Follow the directions as outlined in the user prompt. Then click OK when ready.
38. Once all plates are removed from the deck and sealed, place into a clipseal plastic bag. Click "**OK**" to proceed to the Amphyl wash step.
Note: Before placing the supernatant storage plate into a clipseal bag, check the plate for the transfer of beads. To do this, place the plate onto the ABI magnet (if beads are present they will settle to the bottom of the plate) and then inspect the bottom of the plate. If beads are present refer to the Section 15, Troubleshooting.
39. A final message will advise that the run has completed. Click "**OK**".

9.8. Finalising the MP II Run

1. Transfer left over Resin solution from the 12 channel plate and the Lysis Buffer-DTT (wearing safety glasses) into the glass Lysis-DTT bottle previously used. Discard the 12 channel plate in the biohazard waste bin. Take the bottle to room 6122 and transfer left over reagents into the brown Winchester bottle located in the fume hood.
2. Discard the contents of the Amphyl wash station into the sink and rinse the troughs with Nanopure water.
3. Remove all labware from the deck and clean with 5% TriGene™ followed by 70% Ethanol, and setup for the next run if necessary. In addition, clean the work area around the MP11 instrument with 10% (v/v) Bleach 7x and 70% Ethanol.
4. Remove the tip chute and funnel, rinse with warm tap water to remove any residue inside the chute before cleaning with 5% TriGene and 70% Ethanol.
5. Move the platemap to **C:\PACKARD\EXT PLATE MAPS\Completed Extractions**.

9.9. Importing MP II Log File into AUSLAB

1. Click on the Microsoft Access icon in the WinPrep® main menu to open the MultiPROBE log database.
2. Click on the relevant run Test ID in the Test Selection box. In the Report/Query/Action Selection dropdown menu, select "Report: Test Summary (Sorted by Destination Rack ID)"
3. In the Output Selection dropdown menu, select "File". Save the output file as *.csv format to **C:\PACKARD\EXT PLATE MAPS\EXT LOGS** with the same name as the AUSLAB batch ID and click "**Apply**".
4. Open the log file and check for any errors that may have arisen during the extraction process. Compare the listed errors to any that were encountered during the run. Report any critical errors to the line manager.
5. Copy the log file to **I:\EXTRACTION\EXT A MP11\Logs** or **I:\EXTRACTION\EXT B MP11\Logs** for uploading to AUSLAB.

Automated DNA IQ™ Method of Extracting DNA

6. Import the log file, entering the path, filename and extension (e.g. I:\EXTRACTION\EXT A MPII\Logs\CWIIQEXT20071115_01.csv) and press **[Enter]**. For more details on importing files into AUSLAB refer to Batch functionality in AUSLAB SOP (QIS 24469).

9.10. Importing Extraction “results” into AUSLAB

1. Import the results file, entering the filename and extension. For more details on importing files into AUSLAB refer to Batch functionality in AUSLAB SOP (QIS [24469](#)).
2. The file will be imported into AUSLAB and appear in the DNA file table.
3. Highlight entry and press **[Enter]**, for access to the DNA results table.
4. Page down through the table and check that all sample results have been imported.
5. Press **[SF8] Table Sort Order**, this sorts the table, sorting samples that have failed Autovalidation to the top. Samples that have failed are those that have processing comments present.
6. For samples that have failed, check the **Processing Comments**, by entering into the sample.
7. If processing comments state sample is to be sent to another batch type **other** than quant. Proceed with the following steps:
 - a. Request the appropriate rework test code via the **[SF7]** results history table and the **[SF8]** request rework functions (e.g. samples requiring Microcon, NucleoSpin and pooling).
 - b. Press **[Esc]** to exit back to the DNA results table. Do not toggle accept.
 - c. Add the extraction batch ID into the 9PLEX or 9FTAR completed date fields in AUSLAB.
8. If processing comments do not state next step the sample will be processed as normal:
 - a. Press **[Esc]** to exit back to the DNA results table.
 - b. Highlight any entries to be changed and press **[SF7] Toggle Accept**.
9. Press **[F7] Complete Batch**, all of the samples with a **Yes** in the Accept column will be transferred to the outstanding Quant Casework or Quant Reference lists.
10. File the extraction worksheet into the relevant folder in Room 6117.

9.11. Sample Storage

Refer to “Analytical Sample Storage” (QIS 24255) for how to store the DNA extract Nunc™ Bank-It™ tubes, ABgene 96-deep well and Axygen store plates.

10. TROUBLESHOOTING

1. If the resin is not pipette mixing correctly (eg. resin is not being drawn up into the pipette tip), pause the Winprep program and check that the support tile is clicked onto the MPII deck correctly. Alternatively, pipette and mix resin manually one more time from the corresponding column in the 12 channel plate.
2. If the Gripper is not placing the rack/plate properly on the deck, pause the program and manually adjust the rack/plate properly on the adapter support tile.
3. If the Gripper appears to be not gripping a plate correctly, or the pipetting alignment of a labware appears to be incorrect, the deck and / or labware can be re-calibrated. This is covered in the Operation and Maintenance of the MultiPROBE® II PLUS HT EX and MultiPROBE® II PLUS HT EX with Gripper™ Integration Platform SOP (QIS [23939](#))

11. VALIDATION

- Nurthen, T., Hlinka, V., Muharam, I., Gallagher, B., Lundie, G., Iannuzzi, C. “Project 9. Automation: Report on the Evaluation of DNA Extraction Chemistries.” June 2007.

Automated DNA IQ™ Method of Extracting DNA

- Nurthen, T., Hlinka, V., Muharam, I., Gallagher, B., Lundie, G., Iannuzzi, C. "Project 11: Report on the Validation of the Automated Extraction Chemistry Kit using the MultiPROBE® II PLUS HT EX with Gripper™ Integration Casework Platform." 2007.
- Nurthen, T., Hlinka, V., Muharam, I., Gallagher, B., Lundie, G., Iannuzzi, C. "Project 13: Report on the Verification of the Automated Extraction Chemistry Kit using the MultiPROBE® II PLUS HT EX with Gripper™ Integration Casework Platform." 2007.
- Muharam I., Hlinka V., Gallagher, B., Lundie, G., Iannuzzi, C., Nurthen T., McNevin A., Ientile V. "Project 21: A Modified DNA IQ™ Method Consisting of Off-Deck Lysis to Allow Supernatant Retention for Presumptive Identification of α -Amylase" 2008
- Muharam I., Hlinka V., Gallagher, B., Lundie, G., Iannuzzi, C., Nurthen T., McNevin A., Ientile V. "Project 22: A Modified DNA IQ™ Method for Off-Deck Lysis Prior to Performing Automated DNA Extraction" 2008

12. QUALITY ASSURANCE/ACCEPTANCE CRITERIA

- Positive and negative (reagent blank) controls are included in each extraction batch as per Table 6. These controls are processed as normal samples through to completion.
- If any results are obtained from the reagent blank, either at quantification or during the CE QC check, then the possible source of the contamination is investigated. The samples extracted with this control are thoroughly checked and repeated if possible.

13. REFERENCES

1. Boom, R., Sol, C.J.A., Salimans, M.M.M., Jansen, C.L., Wertheim-van Dillen, P.M.E., & van der Noordaa, J., Rapid and Simple Method for Purification of Nucleic Acids. *J. Clin. Microbiol.*, 1990. 28: p. 495-503.
2. Chen., C.W.T.J., C.A., Recovery of DNA Segments from Agarose Gels. *Anal Biochem.*, 1980. 101: p. 339-341.
3. Cowan, C., The DNA IQ™ System on the Tecan Freedom EVO® 100 Profiles in DNA. *Profiles in DNA*, 2006: p. 8-10.
4. Eminovic, I., Karamehić, J., Gavrankapetanović, F. & Heljić, B., A Simple method of DNA Extraction in Solving Difficult Criminal Cases. *MEDARH*, 2005. 59(1): p. 57-58.
5. Greenspoon, S.B., J., Robotic Extraction of mock sexual assault samples using the Biomek® 2000 and the DNA IQ™ System. *Profiles in DNA*, 2002. 5: p. 3-5.
6. Komonski, D.I., Marignani, A., Richard, M.L., Frappier, J.R.H., & Newman, J.C., Validation of the DNA IQ™ System for use in the DNA extraction of high volume forensic casework. *Can.Soc.Forensic Sci.J.*, 2004. 37(2): p. 103-109.
7. Mandrekar, P., V., Flanagan, L., & Tereba, A., Forensic Extraction and Isolation of DNA Form Hair, Tissue and Bone. *Profiles in DNA*, 2002: p. 11.
8. Mandrekar, P.V., Kreneke, B. E., & Tereba, A., DNA IQ™: The Intelligent Way to Purify DNA. *Profiles in DNA*, 2001: p. 16.
9. Marko, M.A., Chipperfield, R., & Birnboim, H.C., A Procedure for the Large Scale Isolation of Highly purified Plasmid DNA using alkaline extraction and binding to glass powder. *Anal. Biochem.*, 1982. 121: p. 382-387.
10. Melzak, K.A., Sherwood, C.S., Turner, R.F.B. & Haynest, C.A., Driving forces for DNA Adsorption to Silica in Perchlorate Solutions. *J. Colloid. Interface Sci.*, 1996. 181: p. 635-644.
11. PerkinElmer, Automated DNA IQ™ System for Mixed Casework Sample DNA Isolation. *MultiPROBE II Liquid Handling - Forensic Workstation Application Guide*, 2004: p. 1-25.
12. Promega, FAQs –DNA IQ™ System.
13. Promega, *Protocols & Applications Guide*. Chapter 9. rev. 7/06.
14. Promega, DNA IQ™ System -Small Casework Protocol. *Promega Technical Bulletin #TB296* 2006. Rev 4/06: p. 1-14.
15. Promega, DNA IQ™ System-Database Protocol. *Promega Technical Bulletin #TB297*, 2006. Rev 4/06: p. 1-14.

Automated DNA IQ™ Method of Extracting DNA

16. Promega, Tissue and Hair Extraction Kit (for use with DNA IQ™) Protocol. Promega Technical Bulletin #TB307, 2006. Rev 5/06: p. 1-11.
17. Promega Corporation 2006 Material Safety Data Sheet. Lysis Buffer. Article number: A826.
18. Schiffner, L.A., Bajda, E. J., Prinz, M., Sebestyen, J., Shaler, R. & Caragine, T.A., Optimisation of a Simple, Automatable Extraction Method to Recover Sufficient DNA from Low Copy Number DNA Samples for Generation of Short Tandem Repeat Profiles. Croat Med J, 2005. 46(4): p. 578 -586.
19. Vandenberg, N., van Oorchot., R. A. H., & Mitchell, R. J., An evaluation of selected DNA extraction strategies for short tandem repeat typing. Electrophoresis, 1997. 18: p. 1624-1626.

14. STORAGE OF DOCUMENTS

All worksheets are stored in the Analytical area (Room 6117).

15. ASSOCIATED DOCUMENTS

- QIS 17120 Operational Practices in the DNA Dedicated Laboratories
 QIS [17171](#) Method for Chelex Extraction
 QIS 17165 Receipt, Storage and Preparation of Chemicals, Reagents and Test Kits
 QIS [23939](#) Operation and Maintenance of the MultiPROBE® II PLUS HT EX and MultiPROBE® II PLUS HT EX with Gripper™ Integration Platform
 QIS 24255 Analytical Sample Storage
 QIS [24256](#) Sequence Checking with the STORstar Instrument
 QIS 24469 Batch functionality in AUSLAB
 QIS [24919](#) DNA Analysis Workflow Procedure

16. AMENDMENT HISTORY

Version	Date	Author/s	Amendments
R0	23 Oct 2007	B. Gallagher, T. Nurthen, C. Iannuzzi, V. Hlinka, G. Lundie, I Muharam.	First Issue
R1	12 Dec 2007	M Harvey, C Iannuzzi, A McNevin	Reviewed and updated after initial training
R2	19 March 2008	M Harvey, B Andersen, C Iannuzzi, A McNevin	Addition of Off-deck Lysis procedure, Retention of fully automated method as Appendix, addition of reagent record tables into Appendix
R3	April 2008	QIS2 Migration Project	Headers and Footers changed to new CaSS format. Amended Business references from QHSS to FSS, QHPSS to CaSS and QHPS to Pathology Queensland
4	13 March 2009	QIS2 migration	Version incremented by one on migration to QIS2
5	03 June 2009	M Aguilera, B Micic, C Iannuzzi, A. Cheng, V. Hlinka, I. Muharam, G. Lundie, C. Weber	Major changes to reflect new procedure. Updated to reflect changes in procedure as an outcome of internal and external audits. Created ver.6.4 ODL in MPII Platforms. Minor changes in procedures using 4titude 4seal

Automated DNA IQ™ Method of Extracting DNA

			heat sealer to seal plates.
6	29 June 2009	A McNevin, K Lancaster	Removed references to retaining lysate and beads, fixed minor formatting errors. Created ver6.5 ODL in MPII Platforms. Substrates now to be retained in 2mL tube

Not Current

17. APPENDIX

17.1. Manual method for extraction using DNA IQ™

17.1.1. Sampling and Sample Preparation

Refer to section 9 above.

17.1.2. QC samples

All extraction batches require two controls to be registered. The registration of control samples is covered in the DNA Analysis workflow procedure (QIS [24919](#))

17.1.3. Creating the Extraction Batch and Locating Samples

Refer to “DNA Analysis Workflow Procedure” (QIS 24919).

17.1.4. Procedure (No Retain Supernatant)

1. Turn on the Eppendorf Thermomixer and set the temperature to 37°C. Alternatively an appropriately calibrated hot block may be used.
2. Remove DTT 1M (small aliquot) and Prot K 20mg/mL (small aliquot) from freezer to thaw.
3. Label for each sample:
 - Original sample tube; 2mL SSI tube (if original sample is not in a 2mL tube); Spin basket or 2mL tube; and Nunc™ Bank-It™ storage tube.

Note: Spin baskets are not required for the Negative Extraction control, tape lifts, nails and other non absorbent substrates. For these samples, excluding the Negative Extraction control, label a 2mL tube instead of a spin basket. Substrates will be retained into original 1.5mL or 2mL after being processed in a spin basket.
4. Have a second operator perform a sequence check of all tubes. This person must enter the sequence check in AUSLAB.
5. Using Table 7, prepare Extraction Buffer, Lysis Buffer & Resin solution. Ensure that the DNA IQ™ Resin solution is thoroughly vortexed prior to use.

Note: Reagents need to be prepared fresh before each run and Lysis Buffer-DTT solution and Resin solution need to be prepared in the fume hood.

Table 7. Table of reagent volumes for DNA IQ Manual Extraction

Reagent (volume per sample)	Constituent	Volume per sample (µL)	Volume for 12 Samples (mL)	Volume for 24 Samples (mL)
Extraction Buffer (300µL/sample)	TNE Buffer	277.5	4.0	8.0
	Prot K (20mg/mL)	15.0	0.216	0.432
	Sarcosyl (40% w/v)	7.5	0.108	0.216
Lysis Buffer – DTT (726µL/sample)	Lysis Buffer	660	10.0	20.0
	DTT	6.6	0.1	0.2
Resin-Lysis Solution (50µL/sample)	Lysis Buffer with DTT (from above)	43	0.645	1.29
	DNA IQ RESIN	7	0.105	0.210
DNA IQ 1X Wash Buffer (300µL/sample)	See Reagent preparation		4.0	8.0
DNA IQ Elution Buffer (100µL/sample)	Use directly from Kit		1.4	2.8

Note: The volume of Lysis-DTT Buffer calculated includes the volume used in the Resin-Lysis solution preparation.

Automated DNA IQ™ Method of Extracting DNA

6. Add 300µL of Extraction Buffer. Ensure that large substrates including tape lifts are fully submerged. Vortex, then incubate at 37°C on the Thermomixer at 1000 rpm for 45 minutes. If no thermomixer available, incubate samples within a hotblock, vortex mix at the start and end of the incubation and at least one during the incubation.
7. Remove from the Thermomixer and increase the temperature on the Thermomixer to 65°C (for use in the Elution steps).
8. Transfer substrate to spin basket if required or transfer the lysate to an appropriately labelled 2.0mL tube. Retain original tube containing the substrate if no spin basket used.
9. Centrifuge spin basket at maximum speed (14,000rpm) for 2 minutes.
10. Retain the spin basket and transfer the flow through back into sample tube. Transfer the substrate into a labelled 2mL tube.
Note: If original sample tube is not a 2mL tube, transfer flow through from spin basket and the supernatant from the original tube into a 2mL tube.
11. Add 550µL of Lysis-DTT Buffer solution.
12. Add 50µL of DNA IQ™ Resin-Lysis solution. Pipette mix the Resin-Lysis solution between each addition to keep the resin suspended in solution.
13. Vortex and place on the Multitube shaker set at 1200 rpm for 5 minutes.
14. Remove from the Multitube shaker, vortex and immediately place into the magnetic stand. Separation will occur instantly.
Note: If resin does not form a distinct pellet on the side of the tube, or if the magnetic beads have re-suspended while in the stand, vortex the tube and quickly place back in the stand.
15. Carefully remove and discard the solution, ensuring that the resin is not disturbed. Remove from the magnetic stand.
Note: If some resin is drawn up in the pipette tip, gently expel resin back into tube to allow re-separation.
16. Add 125µL of Lysis-DTT Buffer solution. Vortex and place into the magnetic stand.
17. Once separation has occurred, remove and discard the Lysis-DTT Buffer. Remove from the magnetic stand.
18. Add 100µL of 1X Wash Buffer, vortex and place into the magnetic stand. Once separation has occurred remove and discard the Wash Buffer. Remove from the magnetic stand.
19. Repeat the Wash Buffer step (step 18) two times for a total of three washes. Ensure that all of the solution has been removed after the last wash. Remove from the magnetic stand.
20. In a Biohazard hood, uncap the tubes and place the lids down onto a clean rediwipe. Air-dry the resin in the hood for 15 minutes and then recap tubes and remove from the Biohazard hood.
Note: Do not dry for more than 20 minutes, as this may inhibit the elution of DNA.
21. Add 50µl of Elution Buffer by carefully pipetting the solution to the side of the tube, above the pellet. **Do not mix.**
22. Incubate in the Thermomixer at 65°C for 3 minutes and then continue to incubate for a further 3 minutes shaking at 1100 rpm. If no thermomixer available, incubate samples

Automated DNA IQ™ Method of Extracting DNA

within a hotblock, vortex mix at the start and end of 2nd 3 minute incubation. Remove samples.

23. Vortex and immediately place into the magnetic stand while hot to ensure maximum DNA yield during elution.
24. Carefully transfer the DNA extract to the corresponding labelled Nunc™ Bank-It™ tube.
25. Remove from the magnetic stand and repeat the Elution Buffer steps (step 21-24). The final volume after the double elution is approximately 95µL of DNA extract.
26. DNA extracts are stored in temporary storage in freezer 6117-2 (-20°C) located in the workflow area. Tubes containing the original substrate are to be stored in spin basket boxes in freezer 6117-5 located in the workflow area.

17.1.5. Procedure (Retain Supernatant)

1. Turn on the Eppendorf Thermomixer and set the temperature to 37°C. Alternatively an appropriately calibrated hot block may be used.
2. Remove DTT 1M (small aliquot) and Prot K 20mg/mL (small aliquot) from freezer to thaw.
3. Label for each sample: Original sample tube; 2mL SSI tube (if original sample is not in a 2mL tube); 1.5mL tube (for supernatant) these tubes should not be in contact with the substrate; Spin basket or 2mL tube; an extra 2mL tube for spin baskets; Nunc™ Bank-It™ storage tube.
Note: Spin baskets are not required for the Negative Extraction control, tape lifts, nails and other non absorbent substrates. For these samples, excluding the Negative Extraction control, label a 2mL tube instead of a spin basket. Substrates will be retained into original 1.5mL or 2mL after being processed in a spin basket.
4. Have a second operator perform a sequence check of all tubes. This person must enter the sequence check in AUSLAB.
5. Using Table 8, prepare Lysis Buffer & Resin solution. Ensure that the DNA IQ™ Resin solution is thoroughly vortexed prior to use.
Note: Reagents need to be prepared fresh before each run and Lysis Buffer-DTT solution and Resin solution need to be prepared in the fume hood.

Table 8. Table of reagent volumes for DNA IQ Manual Extraction

Reagent (volume per sample)	Constituent	Volume per sample (µL)	Volume for 12 Samples (mL)	Volume for 24 Samples (mL)
Lysis Buffer – DTT (726µL/sample)	Lysis Buffer	660	10.0	20.0
	DTT	6.6	0.1	0.2
Resin-Lysis Solution (50µL/sample)	Lysis Buffer with DTT (from above)	43	0.645	1.29
	DNA IQ RESIN	7	0.105	0.210
DNA IQ 1X Wash Buffer (300µL/sample)	See Reagent preparation		4.0	8.0
DNA IQ Elution Buffer (100µL/sample)	Use directly from Kit		1.4	2.8

Note: The volume of Lysis-DTT Buffer calculated includes the volume used in the Resin-Lysis solution preparation.

6. Add 450µL of TNE buffer and vortex.

Automated DNA IQ™ Method of Extracting DNA

7. Incubate at room temperature for 30 minutes.
8. Vortex, then centrifuge at maximum speed (14,000rpm) for 3 minutes.
9. Remove 150µL of supernatant and place into the respective 1.5mL tube labelled with "sup" (for further testing).
10. Add 14µL of 20mg/mL Proteinase K and 7µL 40% (w/v) Sarcosyl to the original sample tube containing substrate and TNE buffer, then vortex.
11. Vortex, then incubate at 37°C on the Thermomixer at 1000 rpm for 45 minutes. If no thermomixer available, incubate samples within a hotblock, vortex mix at the start and end of the incubation and at least one during the incubation.
12. Remove from the Thermomixer and increase the temperature on the Thermomixer to 65°C (for use in the Elution steps).
13. Transfer substrate to spin basket if required or transfer the lysate to an appropriately labelled 2.0mL tube. Retain original tube containing the substrate in if no spin basket used.
14. Centrifuge spin basket at maximum speed (14,000rpm) for 2 minutes.
15. Retain the spin basket and transfer the flow through back into sample tube.
Note: If original sample tube is not a 2mL tube, transfer flow through from spin basket and the supernatant from the original tube into a 2mL tube. Transfer the substrate into a labelled 2mL tube.
16. Add 550µL of Lysis-DTT Buffer solution.
17. Add 50µL of DNA IQ™ Resin-Lysis solution. Pipette mix the Resin-Lysis solution between each addition to keep the resin suspended in solution.
18. Vortex and place on the Multitube shaker set at 1200 rpm for 5 minutes.
19. Remove from the Multitube shaker, vortex and immediately place into the magnetic stand. Separation will occur instantly.
Note: If resin does not form a distinct pellet on the side of the tube, or if the magnetic beads have re-suspended while in the stand, vortex the tube and quickly place back in the stand.
20. Carefully remove and discard the solution, ensuring that the resin is not disturbed. Remove from the magnetic stand.
Note: If some resin is drawn up in the pipette tip, gently expel resin back into tube to allow re-separation.
21. Add 125µL of Lysis-DTT Buffer solution. Vortex and place into the magnetic stand.
22. Once separation has occurred, remove and discard the Lysis-DTT Buffer. Remove from the magnetic stand.
23. Add 100µL of 1X Wash Buffer, vortex and place into the magnetic stand. Once separation has occurred remove and discard the Wash Buffer. Remove from the magnetic stand.

Automated DNA IQ™ Method of Extracting DNA

24. Repeat the Wash Buffer step (step 23) two times for a total of three washes. Ensure that all of the solution has been removed after the last wash. Remove from the magnetic stand.
25. In a Biohazard hood, uncap the tubes and place the lids down onto a clean rediwipe. Air-dry the resin in the hood for 15 minutes and then recap tubes and remove from the Biohazard hood.
Note: Do not dry for more than 20 minutes, as this may inhibit the elution of DNA.
26. Add 50µl of Elution Buffer by carefully pipetting the solution to the side of the tube, above the pellet. **Do not mix.**
27. Incubate in the Thermomixer at 65°C for 3 minutes and then continue to incubate for a further 3 minutes shaking at 1100 rpm. If no thermomixer available, incubate samples within a hotblock, vortex mix at the start and end of 2nd 3 minute incubation. Remove samples.
28. Vortex and immediately place into the magnetic stand while hot to ensure maximum DNA yield during elution.
29. Carefully transfer the DNA extract to the corresponding labelled Nunc™ Bank-It™ tube.
30. Remove from the magnetic stand and repeat the Elution Buffer steps (step 26-29). The final volume after the double elution is approximately 95µL of DNA extract.
31. DNA extracts & retained supernatants (“sup” tubes) are stored in temporary storage in freezer 6117-2 (-20°C) located in the workflow area. Tubes containing the original substrate are to be stored in spin basket boxes in freezer 6117-5 located in the workflow area.

17.1.6. Sample storage

Refer to “DNA Analysis Workflow Procedure” (QIS [24919](#)).